

Course Instructor: Md. Nazmul Abedin Khan
Course Code: 109
Course Title: CHE


Student name: B M Shahria Alam
ID: 2021-3-60-016
Section: 5

Experiment-5: Determination of the Molar Mass of an unknown monoprotic acid, HA (0.60g in 100mL).

Theory:

A monoprotic acid is an acid that donates only one proton or hydrogen atom per molecule to an aqueous solution, such as: Hydrochloric acid (HCl), Nitric acid (HNO₃), Acetic acid (CH₃COOH). The molar mass of a chemical compound is the mass of a sample of that compound divided by the amount of substance in that sample, measured in moles. Molar masses are usually expressed in g/mol.

Chemical reactions between acids and bases are used to analyze the quantity of a pure substance in a mixture. In this experiment, the standardized base, NaOH, is used to determine the molar mass of an unknown weak acid. Based on the molar mass, you will determine the identity of the unknown weak acid.



Since 1 mol of NaOH reacts with 1 mol of HA, the following expression can be written:

$$(M_{\text{HA}} \times V_{\text{HA}}) = (M_{\text{NaOH}} \times V_{\text{NaOH}}) \dots\dots\dots(1)$$

Where,

M_{NaOH} = Molarity of NaOH

V_{NaOH} = Volume of NaOH

V_{HA} = Volume of HA

M_{HA} = Molarity of HA

Procedure:

a) Dilution of HA:

Retrieve the **1M CH₃COOH** solution and a 100mL volumetric flask from the stockroom. In order to make **0.1M CH₃COOH** solution, take 10mL of **1M** stock **CH₃COOH** solution, transfer it to 100mL volumetric flask and dilute it up to the mark with water.

[Using the formula $M_1 \times V_1 = M_2 \times V_2$; where “1” represents the concentrated conditions (i.e. stock solution Molarity, 1M and volume) and “2” represents the diluted conditions (i.e. desired volume, 100mL and Molarity, 0.1M)].

b) Determination of the concentration of HA:

Take 10 mL of HA (CH₃COOH, as 0.60g in 100 mL) in a conical flask with the help of pipette. Add two drops (0.1mL) of phenolphthalein indicator. Titrate until the last drop of NaOH solution leaves a permanent pink color in the solution and record the final reading in Table.

Calculate the difference between two burette readings (initial and final), which is the amount of

NaOH required neutralizing CH₃COOH. [always subtract volume of CH₃COOH (10mL) and volume of indicator (0.1mL) from total volume, which was added before titration]

Data:

Volume of NaOH, ml	pH of Acid
0	2.88
0.2	3.19
0.4	3.42
0.6	3.58
0.8	3.71
1	3.81
1.2	3.9
1.4	3.97
1.6	4.04
1.8	4.1
2	4.16
2.2	4.21
2.4	4.26
2.6	4.3
2.8	4.35
3	4.39
3.2	4.43
3.4	4.47
3.6	4.51
3.8	4.55
4	4.58
4.2	4.62
4.4	4.65
4.6	4.69
4.8	4.72
5	4.76
5.2	4.79
5.4	4.83
5.6	4.86
5.8	4.9
6	4.93
6.2	4.97
6.4	5.01
6.6	5.05
6.8	5.08
7	5.13
7.2	5.17
7.4	5.21
7.6	5.26
7.8	5.31
8	5.36
8.2	5.42

8.4	5.48
8.6	5.55
8.8	5.62
9	5.71
9.2	5.82
9.4	5.95
9.6	6.14
9.8	6.45
10	8.7
10.2	10.98
10.4	11.28
10.6	11.45
10.8	11.58
11	11.67
11.2	11.74

Calculation:


Va x Ma = Vb x Mb	
Va, Volume of CH ₃ COOH	10 ml
Ma, Molarity of CH ₃ COOH	?
Vb, Volume of NaOH	10.2 ml
Mb, Molarity of NaOH	0.1M
Ma=Vb*Mb/Va	0.102M

Percentage of error
Error = (Theo.value-Exp.value)/(Theo.value) ×100%
= [(0.1-0.102)/0.1]*100%
= 2%

(0.6g in 100ml)		
1000ml	1M	1 mol
100 ml	1M	100/1000
100 ml	0.102M	(100*0.102)/1000
		0.0102 mole

M	= m/n
	= 0.6/0.012
	= 59g/mole

Titration curve



VIRTUAL LAB: Strong Acid and Base Problems

We are pleased to announce a new HTML5 based version of the virtual lab. Please use Firefox or Chrome web browser to access this page, errors have been reported when using Internet Explorer.

[Introductory Video and Support Information](#)

Virtual Lab
File ▾ Edit ▾ View ▾ Help ▾
EN Strong Acid Textbook Problems

Stockroom +

Information ≡

H ⁺	1.80527e-12
OH ⁻	0.00559281
CH ₃ COOH	4.82047e-9
CH ₃ COO ⁻	0.0467290
PhenolphthaleinH	2.85693e-7
Phenolphthalein ⁻	0.0000146676
Na ⁺	0.0523364

Display Absorbance

Download Absorbance Table

Temperature: 25.00°C

25.0 deg

pH: 11.74

Workbench 1

1M CH₃COOH
90.000 mL @ 25.0°C

0.1M CH₃COOH
90.000 mL @ 25.0°C

0.1M NaOH
50.000 mL @ 25.0°C

Phenolphthalein
99.800 mL @ 25.0°C

Unknown CH₃COOH
21.400 mL @ 25.0°C

50 mL Burette
38.800 mL @ 25.0°C

