

Course Instructor: Md. Nazmul Abedin Khan
Course Code: 109
Course Title: CHE

Student name: B M Shahria Alam
ID: 2021-3-60-016
Section: 5


Experiment-4: Standardization of a strong acid (HCl) with a standard weak base (Na₂CO₃).

Theory:

Sodium carbonate is the salt of a weak acid, carbonic acid (H₂CO₃). It is a primary standard substance. Its reaction with HCl is a two-step process as follows:



The overall reaction:



Since 1 mol of Na₂CO₃ reacts with 2 mol of HCl, the following expression can be written:

$$V_a M_a = 2V_b M_b \dots\dots\dots(4)$$

Where,

M_b = Molarity of Na₂CO₃

V_b = Volume of Na₂CO₃ = 10mL

V_a = Volume of HCl

M_a = Molarity of HCl =?

Procedure:

a) Dilution of standard 0.1M Na₂CO₃ from 1M Na₂CO₃ solution:

Retrieve the 1M Na₂CO₃ solution and a 100 mL volumetric flask from the stockroom. Take 10mL of 1M stock Na₂CO₃ solution, transfer it to 100mL volumetric flask and dilute it up to the mark with water.

Dilution:

How to prepare 0.1 M HCl solution from 1.0 M HCl solution

Solution 1: Stock solution; HCl, 1 M

Solution 2: this is the target solution. have to prepare 0.1 M HCl in 100 mL volumetric flask.

Equation is $V_1M_1 = V_2M_2$

Solution 1, mother solution, 1.0 M Na₂CO₃

Solution 2, target solution, 0.1 M Na₂CO₃ in 100 mL

$$V_1 M_1 = V_2 M_2$$

$$V_1 = ?$$

$$M_1 = 1 \text{ M}$$

$$V_2 = 100 \text{ mL}$$

$$M_2 = 0.1 \text{ M}$$

$$V_1 = V_2 \times M_2 / M_1$$

Then we get $V_1 = 10 \text{ mL}$

b) Standardization of HCl:

Take HCl (0.1M) solution in a burette and record the initial reading in Table.

Pipette 10 mL of standard 0.1M Na_2CO_3 solution into a 250 mL conical flask. Add 0.1mL (2 drops) of methyl orange indicator to 10mL Na_2CO_3 solution and the solution's color changes to yellow.

Titrate the solution with the given HCl by continuously adding 0.2 mL (1drop) each time until changing its color from yellow to red. At the end point, record the final reading on burette in table.

Calculate the difference between two burette readings (initial and final), which is the amount of HCl required neutralizing Na_2CO_3

Data:

<u>Volume of HCl, ml</u>	<u>pH of base</u>
<u>0</u>	<u>11.64</u>
<u>0.2</u>	<u>11.54</u>
<u>0.4</u>	<u>11.44</u>
<u>0.6</u>	<u>11.35</u>
<u>0.8</u>	<u>11.27</u>
<u>1</u>	<u>11.19</u>
<u>1.2</u>	<u>11.12</u>
<u>1.4</u>	<u>11.06</u>
<u>1.6</u>	<u>11</u>
<u>1.8</u>	<u>10.95</u>
<u>2</u>	<u>10.89</u>
<u>2.2</u>	<u>10.85</u>
<u>2.4</u>	<u>10.8</u>
<u>2.6</u>	<u>10.76</u>
<u>2.8</u>	<u>10.71</u>
<u>3</u>	<u>10.67</u>
<u>3.2</u>	<u>10.63</u>
<u>3.4</u>	<u>10.6</u>
<u>3.6</u>	<u>10.56</u>
<u>3.8</u>	<u>10.52</u>
<u>4</u>	<u>10.49</u>
<u>4.2</u>	<u>10.45</u>
<u>4.4</u>	<u>10.42</u>
<u>4.6</u>	<u>10.38</u>
<u>4.8</u>	<u>10.35</u>

<u>5</u>	<u>10.31</u>
<u>5.2</u>	<u>10.28</u>
<u>5.4</u>	<u>10.24</u>
<u>5.6</u>	<u>10.21</u>
<u>5.8</u>	<u>10.17</u>
<u>6</u>	<u>10.14</u>
<u>6.2</u>	<u>10.1</u>
<u>6.4</u>	<u>10.07</u>
<u>6.6</u>	<u>10.03</u>
<u>6.8</u>	<u>9.99</u>
<u>7</u>	<u>9.95</u>
<u>7.2</u>	<u>9.91</u>
<u>7.4</u>	<u>9.86</u>
<u>7.6</u>	<u>9.82</u>
<u>7.8</u>	<u>9.77</u>
<u>8</u>	<u>9.71</u>
<u>8.2</u>	<u>9.66</u>
<u>8.4</u>	<u>9.6</u>
<u>8.6</u>	<u>9.53</u>
<u>8.8</u>	<u>9.45</u>
<u>9</u>	<u>9.37</u>
<u>9.2</u>	<u>9.26</u>
<u>9.4</u>	<u>9.13</u>
<u>9.6</u>	<u>8.96</u>
<u>9.8</u>	<u>8.71</u>
<u>10</u>	<u>8.33</u>
<u>10.2</u>	<u>7.95</u>
<u>10.4</u>	<u>7.7</u>
<u>10.6</u>	<u>7.53</u>
<u>10.8</u>	<u>7.4</u>
<u>11</u>	<u>7.3</u>
<u>11.2</u>	<u>7.21</u>
<u>11.4</u>	<u>7.14</u>
<u>11.6</u>	<u>7.07</u>
<u>11.8</u>	<u>7.01</u>
<u>12</u>	<u>6.95</u>
<u>12.2</u>	<u>6.9</u>
<u>12.4</u>	<u>6.85</u>
<u>12.6</u>	<u>6.8</u>
<u>12.8</u>	<u>6.76</u>
<u>13</u>	<u>6.72</u>
<u>13.2</u>	<u>6.68</u>
<u>13.4</u>	<u>6.64</u>
<u>13.6</u>	<u>6.6</u>
<u>13.8</u>	<u>6.56</u>
<u>14</u>	<u>6.53</u>
<u>14.2</u>	<u>6.49</u>
<u>14.4</u>	<u>6.46</u>
<u>14.6</u>	<u>6.42</u>
<u>14.8</u>	<u>6.39</u>

<u>15</u>	<u>6.35</u>
<u>15.2</u>	<u>6.32</u>
<u>15.4</u>	<u>6.28</u>
<u>15.6</u>	<u>6.25</u>
<u>15.8</u>	<u>6.21</u>
<u>16</u>	<u>6.18</u>
<u>16.2</u>	<u>6.14</u>
<u>16.4</u>	<u>6.1</u>
<u>16.6</u>	<u>6.06</u>
<u>16.8</u>	<u>6.02</u>
<u>17</u>	<u>5.98</u>
<u>17.2</u>	<u>5.94</u>
<u>17.4</u>	<u>5.9</u>
<u>17.6</u>	<u>5.85</u>
<u>17.8</u>	<u>5.8</u>
<u>18</u>	<u>5.75</u>
<u>18.2</u>	<u>5.69</u>
<u>18.4</u>	<u>5.63</u>
<u>18.6</u>	<u>5.56</u>
<u>18.8</u>	<u>5.49</u>
<u>19</u>	<u>5.4</u>
<u>19.2</u>	<u>5.29</u>
<u>19.4</u>	<u>5.16</u>
<u>19.6</u>	<u>4.97</u>
<u>19.8</u>	<u>4.67</u>
<u>20</u>	<u>3.91</u>
<u>20.2</u>	<u>3.17</u>
<u>20.4</u>	<u>2.88</u>
<u>20.6</u>	<u>2.71</u>
<u>20.8</u>	<u>2.59</u>
<u>21</u>	<u>2.49</u>

Va x Ma = 2Vb x Mb	
Va, Volume of HCl	20.2
Ma, Molarity of HCl	?
Vb, Volume of Na2CO3	10 ml
Mb, Molarity of Na2CO3	0.1M
Ma=2Vb*Mb/Va	0.099M


Percentage of error:

Error = |(Theo.value-Exp.value)/(Theo.value)|×100%

$$= [(0.1-0.099)/0.1]*100\%$$

$$= 1\%$$

Titration Curve



VIRTUAL LAB: Default Virtual Lab Stockroom

We are pleased to announce a new HTML5 based version of the virtual lab. Please use FireFox or Chrome web browser to access this page, errors have been reported when using Internet Explorer.

[Introductory Video and Support Information](#)

Virtual Lab

File Edit View Help

EN Default Lab Setup

Stockroom +

Information ≡

Name: 10 ml 0.1 M Na₂CO₃

Volume: 31.200 mL

Species (aq)	Molarity
H ⁺	0.00321004
OH ⁻	3.14529e-12
Na ⁺	0.0641026
CO ₃ ⁻²	6.64304e-14
HCO ₃ ⁻	0.00000444258
H ₂ CO ₃	0.0320468
MethylOrangeH	0.00000594505
MethylOrange ⁻	4.65205e-7
Cl ⁻	0.0673077

Temperature: 25.00°C

25.0 deg

pH: 2.49

Workbench 1

Distilled H₂O

2910.0 mL @ 25.0°C

0.1M Na₂CO₃

90.000 mL @ 25.0°C

Methyl Orange

99.800 mL @ 25.0°C

1M Na₂CO₃


90.000 mL @ 25.0°C

50 mL Burette

29.000 mL @ 25.0°C

10 ml 0.1 M Na₂CO₃

31.200 mL @ 25.0°C



Copy

Paste

Format Painter

Clipboard

B

I

U

A

Font

F43

✕

✓

f_x

= 1%

	A	B	C	D	E	
52	10	8.33				
53	10.2	7.95				
54	10.4	7.7				
55	10.6	7.53				
56	10.8	7.4				
57	11	7.3				
58	11.2	7.21				
59	11.4	7.14				
60	11.6	7.07				
61	11.8	7.01				
62	12	6.95				
63	12.2	6.9				
64	12.4	6.85				
65	12.6	6.8				
66	12.8	6.76				
67	13	6.72				
68	13.2	6.68				
69	13.4	6.64				
70	13.6	6.6				
71	13.8	6.56				
72	14	6.53				
73	14.2	6.49				
74	14.4	6.46				
75	14.6	6.42				
76	14.8	6.39				
77	15	6.35				
78	15.2	6.32				
79	15.4	6.28				
80	15.6	6.25				
81	15.8	6.21				
82	16	6.18				
83	16.2	6.14				
84	16.4	6.1				
85	16.6	6.06				
86	16.8	6.02				
87	17	5.98				
88	17.2	5.94				
89	17.4	5.9				
90	17.6	5.85				
91	17.8	5.8				
92	18	5.75				
93	18.2	5.69				
94	18.4	5.63				
95	18.6	5.56				
96	18.8	5.49				
97	19	5.4				
98	19.2	5.29				
99	19.4	5.16				
100	19.6	4.97				
101	19.8	4.67				
102	20	3.91				
103	20.2	3.17				
104	20.4	2.88				
105	20.6	2.71				
106	20.8	2.59				
107	21	2.49				