

Course Instructor: Md. Nazmul Abedin Khan

Course Code: 109

Course Title: CHE

Student name: B M Shahria Alam

ID: 2021-3-60-016

Section: 5


Experiment-3: Standardization of a strong base (NaOH) with a standard weak acid, potassium hydrogen phthalate (KHP)

Theory:



Potassium hydrogen phthalate (KHP) is a monoprotic acidic salt (monopotassium salt of phthalic acid, weak acid) with the formula, $\text{KHC}_8\text{H}_4\text{O}_4$. It is often used as a primary standard for acid-base titrations because it is solid and air-stable, making it easy to weigh accurately and not hygroscopic. KHP dissociates completely in water, giving the potassium cation (K^+) and hydrogen phthalate anion (HP^- or phthalate $^-$).

Equation:

And then as a weak acid hydrogen phthalate reacts reversibly with water to give hydronium (H_3O^+) and phthalate ions.



As KHP is monoprotic, 1 mol of NaOH reacts with 1 mol of KHP according to the following equation:



Therefore, the concentration of standardized NaOH can be determined from reaction (1):

$$(\text{M NaOH} \times \text{V NaOH}) = (\text{M KHP} \times \text{V KHP}) \dots\dots\dots (2)$$

Where,

M_b = Molarity of NaOH

V_a = Volume of KHP

V_b = Volume of NaOH

M_a = Molarity of KHP

Data:

Volume of NaOH, ml	pH of acid
0	2.76
0.1	3.5
0.2	3.8
0.3	3.98
0.4	4.11
0.5	4.22
0.6	4.31

0.7	4.38
0.8	4.45
0.9	4.51
1	4.57
1.1	4.62
1.2	4.66
1.3	4.71
1.4	4.75
1.5	4.79
1.6	4.83
1.7	4.87
1.8	4.91
1.9	4.95
2	4.98
2.1	5.02
2.2	5.05
2.3	5.08
2.4	5.12
2.5	5.15
2.6	5.18
2.7	5.21
2.8	5.25
2.9	5.28
3	5.31
3.1	5.35
3.2	5.38
3.3	5.42
3.4	5.45
3.5	5.49
3.6	5.53
3.7	5.56
3.8	5.6
3.9	5.64
4	5.69
4.1	5.73
4.2	5.78
4.3	5.83
4.4	5.89
4.5	5.95
4.6	6.02
4.7	6.1
4.8	6.19
4.9	6.29
5	6.43
5.1	6.62
5.2	6.94
5.3	10.57
5.4	11.8
5.5	12.09

5.6	12.26
5.7	12.38
5.8	12.47


$V_a \times M_a = V_b \times M_b$	
V_a, Volume of KHP	10 ml
M_a, Molarity of KHP	0.5 M
V_b, Volume of NaOH	5.3 ml
M_b, Molarity of NaOH	??
$M_b = V_a \times M_a / V_b$	0.94 M

Percentage of error:

$$\text{Error} = \left| \frac{\text{Theo.value} - \text{Exp.value}}{\text{Theo.value}} \right| \times 100\%$$

$$= [(1 - 0.94) / 1] \times 100\%$$

$$= 6\%$$



Stockroom



Information



H ⁺	3.37888e-13
OH ⁻	0.0298812
KHP	1.67349e-8
KHP ⁻	0.312500
PhenolphthaleinH	9.08099e-11
Phenolphthalein ⁻	2.49092e-8
Na ⁺	0.342381

Display Absorbance

[Download Absorbance Table](#)

Temperature: 25.00°C

pH: 12.47

Workbench 1



0.500M KHP
15.000 mL @ 25.0°C




Unknown NaOH
50.000 mL @ 25.0°C



Phenolphthalein
99.800 mL @ 25.0°C



50 mL Burette
44.200 mL @ 25.0°C



500 mL Flask
16.000 mL @ 25.0°C

